

### **ULTRA-LOW TEMPERATURE FREEZERS**

Now The Legendary Durability Of The New Brunswick Shaker Is Available For Your Storage Needs Too!









# NBS-86°C Freezers

Entrust the Storage of Your Cultures to the Company You Trust Most for Their Growth

**New Brunswick Scientific's** lab equipment is legendary for design innovation, quality construction and long-term durability. Our -86°C freezers are no exception. Already a leading freezer supplier in Europe, our ultra-low temperature freezers have a 20-year track record of proven performance and are now available worldwide.

#### Nine models in two ranges accommodate a variety of needs:

**Innova models** — utilize Vacuum Insulation Panel (VIP) technology to provide up to 30% more storage capacity, without increasing external dimensions.

**Premium models** — are high-quality, full-featured, yet provide a cost saving alternative by utilizing conventional polyurethane insulation.



#### **ENERGY EFFICIENT**

#### **Tollesbury Technology**

reduces time compressor remains on by 40-50%. Provides quick temperature recovery after door opening, optimizes heat removal, lowers energy bills & extends compressor life. (Innova Models only)

- 1 Outer Door Seals
  - flexible, low-temperature seals prevent cold air loss for lower energy bills and less maintenance
- Quiet Powerful Fan reduces noise and maintenance
- Multiple Compartments with insulated inner doors & lids protect samples and improve energy efficiency

#### SAMPLE SECURITY

- 4 Keyed Locks on Outer Doors & Lids prevent unauthorized access
- 5 User-Defined 4-Digit Password prevents unauthorized changes to alarm and temperature setpoints
- 5 Audible and Visual Alarms indicate high & low temperature conditions, power failure, low battery, blocked filter & fault condition
- 6 Lockable Plate over power & alarm switches prevent accidental shutdown

### Automatic Restart with Non-Volatile Memory

returns setpoints to user-programmed levels after power interrupt

#### **Battery Backup**

maintains temperature settings and activates alarms during power outage

Optional CO<sub>2</sub> & LN<sub>2</sub> Backup Systems protect freezer contents

#### **Remote Contacts**

provide connection to an external, user-supplied, security system



#### OTHER DESIGN ADVANTAGES

#### Innova Models

increase storage by up to 30% over freezers of equal size

**Top-Quality Stainless-Steel Interior** prevents oxidation, scratches or rust; easily cleaned

#### Two-Stage Cascade Refrigeration system uses industrial-grade, commercially-available, hermetically-sealed compressors

#### **Unique Automatic Reset**

protects the microprocessor controller from damage by electrical spikes

#### Internal Voltage Stabilizer

evens out voltage fluctuations, maintaining stable power (60 Hz Innova models only excluding U101). External option on Premium and U101 models


#### Two Pass-Through Ports

allow CO<sub>2</sub>/LN<sub>2</sub> injection or addition of probes

### EASY-TO-USE MEMBRANE KEYPAD, CONTROLS & ALARM SYSTEM

#### FRONT-MOUNTED MICROPROCESSOR CONTROL PANEL

• Ergonomically located • Easy to read and operate • Membrane keypad is easily cleaned



#### **CONVENIENCE**

**Front-Mounted Control Panel** easy to read and reach

#### **Easy-to-Open Door**

provides one-handed operation; cam latch ensures secure door closure

#### **Heated Vent**

prevents vacuum formation, allowing door to be reopened quickly

Low Noise & Low Heat Output allow freezer to be placed directly in the lab

#### **Easy Mobility**

all models fit through standard doors and elevators and include heavy-duty casters. (Compact U101 provided with rear rollers)

**Remote Control & Data Logging** optional CryoCommand™ software with RS-485 interface controls up to 30 freezers

#### **Optional Data Loggers**

from chart recorders to paperless systems simplify validation

#### **Optional Inventory Racking**

available to optimize your storage space

Front-Mounted, Washable Filter easily accessed without tools

#### Stainless-Steel Shelves provided in all uprights.

Adjustable height, high quality

#### **Easily-Removed Inner Doors**

lift off hinges without tools for defrosting & cleaning. (Door hinges in Premium models unscrew)

#### WHAT IS VACUUM INSULATION AND WHY IS IT ADVANTAGEOUS?


Vacuum insulation is an advanced thermal insulation technology that significantly outperforms insulating materials such as closed-cell foams (urethane), foam beads or fiber blankets. These conventional technologies attempt to trap gasses to reduce heat transfer. Vacuum insulation removes the gasses within the insulating space, so the material is placed

under a "vacuum." With air molecules removed, heat transfer is greatly reduced. The result is:

- Thermal insulation performance that far exceeds traditional insulating materials
- Ultra-thin walls for up to 30% more room to store samples

There are other freezers that use Innova Vacuum Insulation Panel technology, but none use NANOGEL™. NANOGEL panels provide:

- The best thermal insulation of any traditional or VIP material on the market
- · An environmentally friendly design, using biodegradable, CFC-free and HCFC-free materials
- A metal-coated polymer film barrier seal to maintain its vacuum. Competitive models are subject to accumulation of gasses over time, resulting in a loss of vacuum level and therefore reduced thermal performance. Even under complete vacuum barrier failure, NANOGEL's performance is two times better than foams.



Powder-Coated Steel Urethane Foam

# INNOVA FREEZERS WITH VACUUM

# UP TO 30% MORE STORAGE CAPACITY — SAME FOOTPRINT

NBS' Innova freezers combine Vacuum Insulation Panel technology with our efficient cooling technology to provide the ultimate in space and power efficiency. These models offer up to 30% more storage capacity than traditionally insulated freezers, without increasing external dimensions. Plus, they come standard with NBS' unique Tollesbury Technology to provide more rapid temperature recovery after door opening for better preserved samples and reduced power consumption.

#### **INNOVA FREEZER ADVANTAGES**

- Greater Internal Storage Capacity in the same footprint through use of Vacuum Insulation Panel technology
- Faster Temperature Recovery after door opening through patented Tollesbury Technology. Improves sample preservation, lowers energy consumption and lengthens compressor life
- Inner Door Seals That Prevent Cold Air Loss and eliminate need for a door heater, for lower energy bills and less maintenance
- And Higher-Grade Seals for outer door and cabinet;
   better-quality internal lids on chest models; and inner doors on uprights that lift off hinges without tools



Keep your samples safe, secure and accessible only by you with Model U101 — a unique, personal-sized upright freezer that fits on or under your lab bench.

#### **CAPACITY & DIMENSIONS**

Model No.	Catalog #	Electrical	Capacity		Max. Box Capacity			Dimensions	Shipping Wt.	
Style	Catalog #	Electrical	L	Cubic Ft.	2"	3″	4"	External	Internal	Shipping wt.
U101 Upright	U9420-0000 U9420-0001	120 V, 60 Hz 230 V, 50 Hz+	101	3.6	60	36	24	32.7" x 35.4" x 22.3" 83 x 90 x 56.6 cm	25.2" x 18.9" x 13" 64 x 48 x 33 cm	256 lbs 116 kg
U535 Upright	U9430-0000 U9430-0002 U9430-0001	120 V, 60 Hz 208/220 V, 60 Hz 230 V, 50 Hz+	535	18.9	336	240	144	76.8" x 31.5" x 34.1" 195 x 80 x 86.7 cm	53.7" x 25.2" x 24.2" 136.5 x 64 x 61.5 cm	551 lbs 250 kg
U725 Upright	U9440-0002 U9440-0001	208/220 V, 60 Hz 230 V, 50 Hz+	725	25.6	504	360	216	76.8" x 40.4" x 34.1" 195 x 102.5 x 86.7 cm	53.7" x 34.1" x 24.2" 136.5 x 86.5 x 61.5 cm	694 lbs 315 kg
C585 Chest	U9400-0000 U9400-0002 U9400-0001	120 V, 60 Hz 208/220 V, 60 Hz 230 V, 50 Hz*	585	20.7	416	288	192	43" x 66.5" x 30.9" 109.2 x 169 x 78.5 cm	30.7" x 47.2" x 24.6" 78 x 120 x 62.5 cm	529 lbs 240 kg
C760 Chest	U9410-0002 U9410-0001	208/220 V, 60 Hz 230 V, 50 Hz*	760	26.9	572	396	264	43" x 80.7" x 30.9" 109.2 x 205 x 78.5 cm	30.7" x 61.4"x 24.6" 78 x 156 x 62.5 cm	628 lbs 285 kg

<sup>(†)</sup> To allow for handles and hinges, add 3.1" (8 cm) to width of upright freezers and 4.3" (11 cm) to the depth of chest freezers.

<sup>(\*)</sup> All 50 Hz freezers are supplied with both U.K. and Schuko cordsets.





#### **INNOVA FREEZER SPECIFICATIONS\***

#### CONSTRUCTION

Insulation: VIP: Ultra-thin Vacuum Insulation Panels with NANOGEL™, encapsulated between a metallic polymer film barrier, bonded to foamed-in-place polyurethane foam; 3.15", 80 mm thick

Exterior: 18 gauge steel,1.2 mm thick. Powder coated, scratch

and rust resistant

Interior: Top-quality, polished 304L stainless steel Heated Air Vent: with plunger prevents vacuum formation Air Filter: Front-mounted; Easily accessed by 1/4 turn screws Access Ports: Two pass-through ports, 3/4", 20 mm dia.

Casters & Feet: Heavy-duty lockable casters with adjustable leveling feet. Two rear rollers provided on Model U101

#### TEMPERATURE CONTROL SYSTEM

Range: Programmable, to -86°C, in 1°C increments, even at ambient temperature up to 32°C

Control: Microprocessor-controlled temperature and alarms with non-volatile memory

Display: Large, bright, digital LED characters, 1/2", 1.3 cm high. Mounted for ergonomic viewing

Optional Remote Control: Factory-installed RS-485 port with optional CryoCommand™ software monitors/controls up to 30 freezers simultaneously

Battery Back-up: Activates alarms and displays temperature during power outages. Separate on-off switch provided for cleaning/servicing/thawing

#### **REFRIGERATION**

Cascade Refrigeration: Hermetically-sealed two-stage cascade system with capacity to cope in high-ambient conditions Compressors: Low-noise, heavy-duty & commercially-available Tollesbury Technology: Patented system reduces temperature recovery time after door opening, minimizing power consumption and optimizing compressor efficiency Circulation: Requires just one fan for reduced power

Refrigerants: CFC-free and HCFC-free

Compressor Lubricants: Biodegradable and commerciallyavailable, high-performance, synthetic polyolester oil.

Prevents "oil-logging"

consumption and noise

#### **SECURITY FEATURES**

Unique 4-Digit Password: For temperature setpoint and alarms, prevents unauthorized changes

Lock and Key: for outer doors/lids. Barrel lock is removable Lockable Security Plate: protects on-off switch and alarm on/off

#### **MONITORING & ALARMS**

S.M.A.R.T. Plus™ Diagnostic Software: Built into the control panel, provides fault codes to trace and solve system errors Remote Alarm Contacts: Provide connection to external usersupplied alarm/monitoring system

Visual Alarms: High/low temperature, power out, system fail, battery low, filter clean and fault analysis. LED lamp for remote control operation

Audible Alarms: High/low temperature, power out and system fail

#### **UPRIGHT MODEL FEATURES**

Door Latch: Injection-molded fiber-reinforced handle with built-in removable lock. Powerful cam-action/reverse camaction latch provides one-handed operation; pulls door firmly closed and pushes door open. Clicks to indicate proper closure Compartments and Shelving: Three compartments with two adjustable-height, stainless-steel shelves. (Model U101: Two compartments, one fixed shelf.) Additional shelving available Insulated Inner Doors: Latchable doors lift off stainless-steel hinges, without tools. Two on U101; three on U535 & U725 Cabinet & Door Seals: Cabinet and outer doors fitted with high-efficiency silicone seals. High-quality seals for inner doors also provided

#### FEATURES UNIQUE TO CHEST MODELS

Door Latch: Positive, single-handed cam action latch with lock Insulated Sub-Lids: Two injection-molded lids on C585; three on C760

#### **CERTIFICATIONS**

**CE Certified and UL Certified** 

(\*) Specifications subject to change without notice (◊) U.S. patent No. 6,449,980. U.K patent No. GB2370874.







# POWER CONSUMPTION, FEATURES & PLUGS

#### **POWER FEATURES & ADVANTAGES**

Internal Voltage Stabilizer (60 Hz Innova Models excluding U101): Evens out high and low voltage fluctuations and optimizes energy usage. Optional external unit available for Premiums and U101

**Automatic Reset:** Protects the microprocessor controller from failure caused by power spikes

**Automatic Restart with Non-Volatile Memory:** Returns alarms and setpoints to user-programmed levels after power interrupt

**Programmed Restart:** Random startup, 1 - 1.5 minutes apart. Prevents power overload

#### PLUG CONFIGURATIONS

Freezers with	Have this	Require this	NEMA Reference			
this Code	Plug	Receptacle	Plug	Receptacle		
А	<b>(1)</b>	0	5-15P	5-15		
В	3	0	5-20P	5-20		
С	⊙	•	6-15P	6-15		
I N	J T E R	NATI	O N A	A L		
D	$\odot$	0	European (Schuko)			
Е	(_)	•	UK Model			

#### POWER CONSUMPTION\*

Brand	Model	Voltage	Amps	Plug Provided	Power Consump. Watts kWh/day		Pull-Down Time /Hrs.
	U101	120 V, 60 Hz 230 V, 50 Hz	10 5	A D/E	530	12.72	3.7
Innova	U535	120 V, 60 Hz ◊ 208/220 V, 60 Hz 230 V, 50 Hz	13 7 7	B C D/E	550	13.20	4.5
IIIIOVa	U725	208/220 V, 60 Hz 230 V, 50 Hz	10 10	C D/E	660	15.84	4.9
	C585	120 V, 60 Hz ◊ 208/220 V, 60 Hz 230 V, 50 Hz	13 7 7	B C D/E	565	13.56	3.5
	C760	208/220 V, 60 Hz 230 V, 50 Hz	10 10	C D/E	665	15.96	5.5
	U410	120 V, 60 Hz <sup>()</sup> 208/220 V, 60 Hz 230 V, 50 Hz	15 7 7	B C D/E	525	12.60	3.75
Premium	U570	120 V, 60 Hz <sup>()</sup> 208/220 V, 60 Hz 230 V, 50 Hz	12 7 7	B C D/E	565	13.56	4.5
	C340	120 V, 60 Hz <sup>()</sup> 230 V, 50 Hz	15 7	B D/E	420	10.1	4.6
	C660	208/220 V, 60 Hz 230 V, 50 Hz	7 7	C D/E	595	14.28	4.0

<sup>(\*)</sup> All models use single phase power. Amperage shown designates "running" amps.

Power Consumption Data & Pull-Down Times are based on an empty freezer set to -80°C, ambient temperature 23°C, using a 240 V, 50 Hz power supply.

#### **INVENTORY RACKING CAPACITIES**

Model	Style	Brand	Cu. Ft.	Compartments per Freezer	Racks per Shelf	Racks per Freezer	Box Size*	Boxes per Rack	Boxes per Freezer	Vials per Freezer**
U101	Upright	Innova	3.6	2	3	6	2" 3"	10 6	60 36	6,000 3,600
U535	Upright	Innova	18.9	3	4	12	2" 3"	28 20	336 240	33,600 24,000
U725	Upright	Innova	25.6	3	6	18	2" 3"	28 20	504 360	50,400 36,000
C585	Chest	Innova	20.7	NA	NA	32†	2″ 3″	13 9	416 288	41,600 28,800
C760	Chest	Innova	26.9	NA	NA	44†	2" 3"	13 9	572 396	57,200 39,600
U410	Upright	Premium	14.5	5	3	15	2" 3"	16 12	240 180	24,000 18,000
U570	Upright	Premium	20.0	5	5	25	2" 3"	16 12	400 300	40,000 30,000
C340	Chest	Premium	12.0	NA	NA	20†	2" 3"	13 9	260 180	26,000 18,000
C660	Chest	Premium	23.3	NA	NA	40†	2" 3"	13 9	520 360	52,000 36,000

<sup>(\*) 4&</sup>quot; Boxes also available. Ask your rep for details

 $<sup>(\</sup>lozenge)$  120 V freezers with this designation require a dedicated 20 Amp line to cover the amps drawn at start-up.


<sup>(\*\*)</sup> Based on a box capacity of 100 vials

<sup>(†)</sup> Half-height racking available

# FREEZER ACCESSORIES Everything you need to preserve the integrity of your samples.



CO<sub>2</sub>/LN<sub>2</sub> controller features selectable setpoint dial. LED lights to indicate when system is activated or valve open



NBS' CryoCommand software



7-day chart recorder



NBS manufactures a variety of racking systems, boxes and dividers to maximize storage space

#### CO<sub>2</sub> & LN<sub>2</sub> CONTROLLER BACK-UP SYSTEM

Safely protect your samples

- Maintain temperature between -55°C and -70°C with Liquid CO<sub>2</sub>, or down to -85°C with LN<sub>2</sub>, subject to environmental conditions
- Coolant addition is pulsed to avoid excessive consumption
- Built-in safety device prevents pulsing while the door is open

## DATA LOGGING and CONTROL CryoCommand™ Software

Remotely log and control temperature and alarms in up to 30 NBS freezers from a single PC

- Easy-to-use screens let you individually change high and low temperature alarms and setpoints
- · Alarm conditions are recorded and signalled by a flashing screen
- Package includes software and factory-installed RS-485 interface

#### Chart Recorder/Data Loggers

Continuously monitor and record temperature for up to 7 days

- Temperature range of +20°C to -90°C, ± 2%
- Front-mounted for easy access. Factory installation recommended
- Contact your sales representative for other options

#### **AUTO-DIALERS**

Auto-dialers can call a set of preprogrammed telephone numbers in the event of an alarm condition and connect right into the freezer's remote alarm port

#### **TEMPERATURE PROBES**

Additional temperature probe for an external alarm system or for validation can be installed upon request

#### **VALIDATION PACKAGES**

Installation and operational qualifications are available. Ask your NBS sales rep for details

#### **INVENTORY RACKING & SHELVES**

Exclusively designed and manufactured for NBS freezers to maximize storage space

- · Anodized aluminum racking is lightweight and durable
- Standard racks hold 2", 3" or 4" boxes
- Custom racking solutions also available
- Waterproof boxes and wax-coated cardboard dividers provide a low-cost alternative to injection-molded boxes
- · Supplemental stainless-steel shelves also available

#### **WARRANTY & AFTER SALES SUPPORT**

#### All NBS freezers are covered by a FIVE-YEAR Limited Warranty:

Two years covering parts and labor — with THREE ADDITIONAL years covering parts. Additionally, VIP insulation panels are warranteed for TWELVE YEARS. Customizable extended warranties and preventive maintenance plans can also be tailored to best suit your individual requirements. Outside of the U.S. and Europe, contact your local NBS sales representative for warranty details.

Every NBS freezer is also equipped with unique S.M.A.R.T. Plus™ diagnostic software to help quickly identify the cause of a fault or setpoint variance. Easily accessed through the front control panel, S.M.A.R.T., Systems Monitoring and Reporting Technology, can be used to identify and correct problems, often over the phone at no charge, with the aid of our trained service technicians. S.M.A.R.T. Plus not only eliminates unnecessary service calls, but provides advance diagnosis, so should an on-site repair become necessary, the service engineer arrives at your facility with the components needed to complete the repair.

Compressor and biodegradable gasses are commercially available, for ready replacement, anywhere in the world.

**NBS** freezers are backed by one of the largest service organizations, comprised of in-house technicians, field engineers, and factory-authorized service centers in your local area.







### **New Brunswick Scientific**

P.O. Box 4005 44 Talmadge Road Edison, NJ 08818-4005 USA

In North America: 800-631-5417 Elsewhere: 732-287-1200 Fax: 732-287-4222

E-Mail: bioinfo@nbsc.com www.nbsc.com

BELGIUM New Brunswick Scientific NV-SA Stationsstraat 180/4 3110 Rotselaar Tel: +32 (0) 16 56 28 31 Fax: +32 (0) 16 57 27 53 E-Mail: sales@nbsnv-sa.be

FRANCE
New Brunswick Scientific S.A.R.L.
3, Rue des Deux - Boules
75001 Paris
Tel: +33 (0) 1-4026 2246
Fax: +33 (0) 1-4026 5423
E-Mail: sales@nbssarl.fr

GERMANY
New Brunswick Scientific GmbH
In der Au 14
D-72622 Nürtingen
Tel: +49 (0) 7022-932490
Fax: +49 (0) 7022-32486
E-Mail: sales@nbsgmbh.de

THE NETHERLANDS
New Brunswick Scientific B.V.
Kerkenbos 1101, 6546 BC
Nijmegen
P.O. Box 6826, 6503 GH Nijmegen
Tel: +31 (0) 24 3717 600
Fax: +31 (0) 24 3717 640
E-Mail: sales@nbsbv.nl

UNITED KINGDOM New Brunswick Scientific (UK) Ltd.

17 Alban Park, Hatfield Road St. Albans, Hertfordshire AL4 0JJ Free Tel: 0800 581331 Tel: +44 (0) 1727 853855 Fax: +44 (0) 1727 835666 E-Mail: bioinfo@nbsuk.co.uk Web: www.nbsuk.co.uk

CHINA, Shanghai New Brunswick Scientific Co., Inc. Room 1501 No. 18 Xiang Jiang Building Lane 1265, Zhongshan Rd. West Shanghai 200051

Tel: +86 21 3223 0203 Fax: +86 21 6278 7182 E-Mail: nbschc@online.sh.cn